

Agilent AFM/Raman System Combines Atomic Force Microscopy and Raman Spectroscopy

Data Sheet

Features and Benefits

- New, fully integrated, high-precision AFM/Raman system
- High performance Raman system for material ID and analysis
- Advanced software enables one-click operation
- Simplified alignment of the Raman laser on the AFM tip
- Tip Enhanced Raman Spectroscopy (TERS) ready system for transparent samples
- Incubator perfusion cell sample plate facilitates dynamic studies in liquids and gases
- Enhanced PicoView software flexibility for plug-ins and scripting

Applications

- AFM/Raman co-localization
- TERS
- AFM force-volume spectroscopy
- Life science and materials science



The Agilent AFM/Raman System.

Overview

This new system seamlessly integrates the Agilent 6000ILM AFM (atomic force microscope) and a HORIBA XploRA INV (inverted Raman microscope), letting researchers go beyond the optical diffraction limit to achieve nanoscale resolution imaging as they perform Raman spectroscopy.

System capabilities include AFM/Raman co-localization, tip-enhanced Raman spectroscopy, and AFM force-volume spectroscopy. The system is ideal for conducting advanced life science research, including studies of cell membranes, the surface structure of cells, single DNA/RNA strands, individual proteins, single molecules, and biopolymers. It is also ideal for investigating novel materials (e.g., graphene).



Agilent Technologies

AFM/Raman co-localization can be performed simultaneously using PicoView to control the Raman Spectroscopy system. Raman mapping is achieved by acquiring a complete Raman spectrum at each pixel of a 2D image to create a detailed chemical image of the sample. This image is generated by plotting the peak intensity (material concentration), peak position (molecular structure or material stress), or peak width (crystallinity). The system precisely overlays the detailed chemical image with a 3D AFM topography image.

An XY piezo stage is used for both Raman and AFM measurements. The 6000ILM AFM provides a wide range of direct sample property measurement (e.g. elasticity) which can be correlated via Raman with chemical composition. Tip-enhanced Raman spectroscopy (TERS) is possible by focusing the Raman laser at the tip of a special AFM probe leading to field enhancement of the laser/sample interaction and thus yielding higher Raman signals as well as increased spatial resolution.

An extensive range of AFM force-volume spectroscopy capabilities is also provided. Results are fast and reliable. Researchers can use their own algorithms and plug-ins, select their own data points on-the-fly, acquire force-curve measurements on any data point, and change experimental parameters in real time.




Figure 1. AFM co-localization and Raman mapping of red blood cells. (A) Raman spectrum. (B) Raman map of peak at 739 cm⁻¹. (C) AFM topography image. (D, E) High resolution AFM topography of the middle cell on the left side in image (B).

Advantages

The 6000ILM AFM utilizes a computer-controlled laser with automated photodetector alignment. A high-stability, low-noise motorized stage directs the movement of the sample beneath the AFM tip for measurement. Easy-to-use PicoView software provides impressive built-in functionality, such as point-and-shoot AFM imaging based on an optical image.

Additional advantages include unrivaled in-liquid AFM imaging via Agilent's patented MAC Mode, the availability of an incubator perfusion cell sample plate to facilitate dynamic studies in liquids and gases, and a top-view video optics package that offers the ability to see opaque samples while scanning.

The XploRA INV from HORIBA combines the exclusive automation features and small footprint of the standard XploRA Raman microscope

with the unique sampling capabilities of an inverted microscope, especially important for demanding biological applications.

PicoView plug-in and scripting flexibility further expand the usefulness of the fully integrated AFM/Raman system.




Figure 2. (A) Raman spectrum: G-band shows visible in graphite and graphene. 2D-band: only visible in graphene. Defect band: visible at the edges of graphene flakes. (B) Raman map: 2D, G and defect band overlaid. (C) AFM: 3D topography.




Figure 3. Simultaneous elasticity and Raman mapping of red blood cells. (A) Optical image of red blood cells; cell marked with red circle used for Raman mapping and elasticity. (B) Raman map. (C) AFM elasticity image.

Agilent 6000 ILM System Specifications

AFM Scanner	
X-Y (Closed Loop) Scan Range	100 μ m x 100 μ m
Z (Closed Loop) Scan Range	30 μ m
Scanner Noise	0.3 nm (X–Y) 0.1 nm (Z)
Scanner Laser	Infrared (980 nm)
Supported Microscopes	
Nikon	Ti-2000 Eclipse Series
Stages	
Stage Travel	4 mm x 6 mm
Motorized Stage Accuracy	\pm 2.5 μ m
System Noise	
	< 0.25 nm RMS (Air)
	< 0.25 nm RMS (Fluid)
Available Sample Plates	
Cover Slip	
Size	22–30 mm dia
Thickness	>170 μ m
Microscope Slide	
Size	1" x 3" (25 mm x 75 mm)
Petri Dish	
Size	35 mm x 10 mm or 50 mm x 9 mm Plastic or plastic with glass bottom

For full 6000 ILM specifications please see 6000 ILM data sheet.

HORIBA XploRA

visit: <http://www.horiba.com/us/en/scientific/products/raman-spectroscopy/raman-spectrometers/raman-microscopes/details/xplo-ra-tm-124/>

AFM Instrumentation from Agilent Technologies

Agilent Technologies offers high-precision, modular AFM solutions for research, industry, and education. Exceptional worldwide support is provided by experienced application scientists and technical service personnel. Agilent's leading-edge R&D laboratories are dedicated to the timely introduction and optimization of innovative and easy-to-use AFM technologies.

www.agilent.com/find/nano

Americas

Canada	(877) 894 4414
Latin America	305 269 7500
United States	(800) 829 4444

Asia Pacific

Australia	1 800 629 485
China	800 810 0189
Hong Kong	800 938 693
India	1 800 112 929
Japan	0120 (421) 345
Korea	080 769 0800
Malaysia	1 800 888 848
Singapore	1 800 375 8100
Taiwan	0800 047 866
Thailand	1 800 226 008

Europe & Middle East

Austria	43 (0) 1 360 277 1571
Belgium	32 (0) 2 404 93 40
Denmark	45 70 13 15 15
Finland	358 (0) 10 855 2100
France	0825 010 700*
	*0.125 €/minute
Germany	49 (0) 7031 464 6333
Ireland	1890 924 204
Israel	972-3-9288-504/544
Italy	39 02 92 60 8484
Netherlands	31 (0) 20 547 2111
Spain	34 (91) 631 3300
Sweden	0200-88 22 55
Switzerland	0800 80 53 53
United Kingdom	44 (0) 118 9276201

Other European Countries:

www.agilent.com/find/contactus

Product specifications and descriptions in this document subject to change without notice.

© Agilent Technologies, Inc. 2013
Printed in USA, May 7, 2013
5991-2338EN



Agilent Technologies