




**WARNING :** When unpacking the analyzer, if the exterior of the analyzer appear to be damaged, **DO NOT turn ON the power switch**, to avoid any electrical shock.


**1** Verify contents of shipment



**Turn ON Line Switch**




**2** Connect the power cable to AC outlet & rear panel of the unit



**4** Press the Standby Switch

The unit is ready for use.  
For detailed instruction,  
read the Installation Guide.



**6** Read the license agreement, then click "Agree" button.



**5** Press <ENTER> to perform Initial Registration of the Windows OS

